

## Dual-Tension Ankle Support

Anatomic shape. Compression class II, closed heel. Circular knitted, thus without seams. Microfiber comfort.


### Indications

For ankle joint injuries, after ankle and ankle joint operations, subluxatio supinatoria, pes valgus troubles, after removal of plaster, loose ligaments, sinew irritations.

### Sizing

Circumference of ankle as shown ( 1cm = 0.3937 in).

Reference	Size	b	Color
GB-0545001	Small	Up to 21 cm	Beige
GB-0545002	Medium	21-23 cm	Beige
GB-0545003	Large	23-25 cm	Beige
GB-0545004	X-Large	25-27 cm	Beige



Measuring point  
*b* = Circumference of ankle corresp. *B* measurement compression stocking

## TaloStabil® Eco Ankle Support

With 2 silicone pads. Anatomical shape. Flat knitted. Compression class II.

Relief of pressure at ankle ending. Extension at side of foot sole avoids curling.

**For left and right application.**


### Indications

After removal of plaster, **distorsions**, arthrosis, achillodynia. Ligament instability.

### Sizing

Circumference of ankle as shown ( 1cm = 0.3937 in).

Reference	Size	Circumference	Color
GB-0546001BE	Small	Up to 21 cm	Beige
GB-0546002BE	Medium	21-23 cm	Beige
GB-0546003BE	Large	23-25 cm	Beige
GB-0546004BE	X-Large	27-27 cm	Beige
GB-0546101BL	Small	Up to 21 cm	Blue
GB-0546102BL	Medium	21-23 cm	Blue
GB-0546103BL	Large	23-25 cm	Blue
GB-0546104BL	X-Large	25-27 cm	Blue



Measuring point  
*b* = Circumference of ankle corresp. *B* measurement compression stocking

## Stabilo®Pren Ankle Support

Neopren bandage with closed heel.  
With 2 Velcro® closures for variable compression and support of supination or pronation.  
Limitation of talus advance.

### Indications

Instability of ankle joint area.



Reference	Size	Shoe size (Men)	Shoe Size (Women)	Color
GB-1822001	X-Small	Up to 6	7½	Black
GB-1822002	Small	6 - 7	7½ - 8	Black
GB-1822003	Medium	7 - 8½	8 - 10	Black
GB-1822004	Large	8½ - 12	10 - 13½	Black
GB-1822005	X-Large	12 - 13½	13½ - 15	Black

## Stabilo® Ankle Support

U-shaped stabilization element made of plastic material and elastic Stabilo® Velcro® strap enable limitation of talus prolapse.  
Comfortable and dry wearing feeling due to naturally breathing and humidity draining material.  
Individual adjustment by lacing.  
Washable.


### Indications

- To stabilize the ankle joint after operations or removal of plaster.
- For immobilization in the case of inflammatory irritations, rheumatism, arthritis, overstrained pes valgus.
- Tendinitides.

### Sizing

Circumference of ankle as shown ( 1cm = 0.3937 in).

Reference	Size	b	Color
GB-1141001	X-Small	Up to 19 cm	Blue/Black
GB-1141002	Small	19 - 21 cm	Blue/Black
GB-1141003	Medium	21 - 23 cm	Blue/Black
GB-1141004	Large	23 - 25 cm	Blue/Black
GB-1141005	X-Large	25 - 27 cm	Blue/Black



Measuring point  
*b* = Circumference of ankle corresp.  
*b* measurement compression stockings

## AchilloStabil® Ankle Support

With especially padded seam. Three-dimensional form-knitted dual-tension elastic bandage. Compression class II. Viscoelastic profile pads at para-achilles area. Pressure relieve of achilles tendon and intermittent massage effect for better blood circulation. Pressure reduced edge. Additional relief of achilles tendon by special heel wedge (optional). With heel cushion to non-affected leg to equalize length differences.


### Indications

Acute and chronic achillodynia, postoperative (follow-up treatment) after ruptures.

### Sizing

Circumference of ankle as shown ( 1cm = 0.3937 in).

Reference	Size	b	Color
GB-0549001	Small	Up to 21 cm	Beige
GB-0549002	Medium	21-23 cm	Beige
GB-0549003	Large	23-25 cm	Beige
GB-0549004	X-Large	25-27 cm	Beige



*Measuring point*

*b = Circumference of ankle corresp. b measurement compression stockings*

## NEURODYN Foot-Lift Orthosis

Fabric Orthosis for the dynamic compensation of flaccid paralysis.

### Product Features:

- Fabric orthosis for dynamic compensation of flaccid paralysis.
- Orthopaedic rehabilitation should compensate for lost functionality without hindering remaining functions.
- Dynamics instead of rigid splints – comfort instead of pressure.
- Free plantar flexion.



Reference (*)	Size	Circumference of leg above ankle	Color (**)
ER-07077-1XX	1	up to 19 CM	Tan
ER-07077-2XX	2	19 - 21 CM	Tan
ER-07077-3XX	3	21 - 23 CM	Tan
ER-07077-4XX	4	23 - 25 CM	Tan
ER-07077-5XX	5	25 - 27 CM	Tan

(\*) To order, replace XX by LT for Left side or RT for Right side

(\*\*) We carry the tan color instead of the black shown on the picture

## Arthrofix Orthosis

More than just an alternative to a plaster cast.

### Product Features:

- This orthosis can be put on without moving the ankle.
- Sliding joint and inner pad for undiminished walking enjoyment.
- Its light construction makes it highly suitable for everyday use.
- Precise fixing due to individual adjustment of the thermoplastic material.



### Indications:

- Conservative therapy for sprained ankles, no matter how severe.
- Post-operative aftercare.

Reference (***)	Size	Shoe Size
ER-07075-SX	1	37 - 39
ER-07075-MX	2	40 - 42
ER-07075-LX	3	43 - 45

(\*\*\*) To order, replace X by L for Left side or R for Right side

## Metatarsal Bandage

Elastic bandage.  
Also available with metatarsal pad  
(No. GB-112070xx below)

**Packing unit:** pair

**Indications**  
For flat foot and splayfoot.

**Sizing**  
Circumference around forefoot as  
shown ( 1cm = 0.3937 in).

Reference	a (*)	Color
GB-11107016	16 cm	Beige
GB-11107017	17 cm	Beige
GB-11107018	18 cm	Beige
GB-11107019	19 cm	Beige
GB-11107020	20 cm	Beige
GB-11107021	21 cm	Beige
GB-11107022	22 cm	Beige
GB-11107023	23 cm	Beige
GB-11107024	24 cm	Beige
GB-11107025	25 cm	Beige
GB-11107026	26 cm	Beige



## Metatarsal Bandage with Met Pad

Elastic bandage. Metatarsal pad inserted in aperture.  
Also available without met pad (No. GB-111070xx  
above)

**Packing unit:** pair

**Indications**  
For flat foot and splayfoot.

**Sizes**  
Circumference around forefoot as

Reference	a (*)	Color
GB-11207016	16 cm	Beige
GB-11207017	17 cm	Beige
GB-11207018	18 cm	Beige
GB-11207019	19 cm	Beige
GB-11207020	20 cm	Beige
GB-11207021	21 cm	Beige
GB-11207022	22 cm	Beige
GB-11207023	23 cm	Beige
GB-11207024	24 cm	Beige
GB-11207025	25 cm	Beige
GB-11207026	26 cm	Beige



(\*) Measuring point: a = Circumference ossa metatarsalia corpus I-V

## Hallux Valgus Brace

Elastic brace with ball protection and splayfoot pad.

**Packing unit:** 1 Item

**Sizes:** Foot circumference (a) from 17 to 26 CM, Left and Right

**Ref.: AF-9046004XXYY (\*)**

(\*) To order, replace in the reference number XX by the foot circumference and YY by LT for the left side, by RT for the right side.

Example:

AF-904600418LT      Left foot, 18CM Circumference



## Hallux Valgus Night Splint

Synthetic clip to correct the ball of the foot and the non-aligned big toe.

**Packing unit:** 1 Item

**Sizes:** Comes in three sizes

Left and Right

Reference(**)	Size	Shoe Size
OD-1333436X	1 - Small	34 - 36
OD-1333740X	2 - Medium	37 - 40
OD-1334146X	3 - Large	41 - 46



(\*\*) To order replace X in the reference with L for left foot or R for right foot

## Valco® Hallux Valgus Bandage

The **anatomically adapted shape** of the splint guarantees the physiological correction of the basic joint of the big toe. **Good wearing comfort** thanks to the generously dimensioned foam paddings which safely protect in the critical area of the toe semi-shell and in the bunion area. To avoid excessive heat, the splint is provided with a **ventilation hole in the toe area**. Two adjustable Velcro® closures allow a precise adjustment of the moments of leverage with the splint remaining in the same **optimum position**.



### Indications

- Hallux valgus,
- Post-operative capsule relief.

### Sizes

Comes in three sizes, for men and women

Reference (*)	Shoe size (US)	Color
GB-930010SX	W 5-8 / M 5	White/blue
GB-930010MX	W 9-11 / M 6-8	White/blue
GB-930010LX	W 12-14 / M 9-13	White/blue

(\*) To order, replace the X in the reference number by R for

## Toe Alignment Splint

Universal toe alignment (one size fits all - left and right), splint. Constructed of soft cotton backed elastic material that will fit into standard footwear. Designed to maintain position of soft tissues until they heal following surgery for bunions or hammer toes. Reduces incidence of recurrence of the deformity following surgical correction.



**Sizes:** One size only.

**Ref.:** ALI-6078      Universal, sold by the unit

## Hammer Toe Correction Bandage

Individually adjustable, light, easy to put on. One bandage is necessary per correction. Several bandages can be worn at the same time.

### Indications

To correct hammer toes.

**Sizes:** One size only.

**Ref.:** GB-112090      Universal, sold by the unit



## Link® Toe Splints

Padded and Thermoplastic

### Indications

Fractures, Distorsion, Mal-Alignment, Hammer Toe, Toe Paralysis, Postoperative Protection, Clubfoot

**Instructions:** After having shaped the splint, fix it with two layers of medical adhesive tape. Avoid applying the tape too tightly in order to avoid edema. The proximal tape wrap is positioned around the ball of the foot and the distal wrap around the affected toe. It is suggested that a sock be worn for additional support of the splint and tape. If necessary, Link® Toe Splints may be reshaped. The material becomes malleable at 65°C (150°F). Use of a hair dryer or submersion into hot water will soften the material enough for further modeling.



Reference	Size	Side
PK-11300001	For Lesser Toe	Right
PK-11300002	For Lesser Toe	Left
PK-11300008	For Big Toe	Right
PK-11300009	For Big Toe	Left
PK-11300005	For Mid Toes	Left and Right



## Silicone Heel Spur Cushion with Softspot

Heel cushion made of silicone with lateral SoftSpot insert. With raised edges. Anatomically designed soft bedding in the heel spur area. Defined and efficient pressure relief by SoftSpot-technology. Effective reduction of shock impact onto ankle, knee and hip joints, as well as onto spine.

**Height:** 6 mm

**Packing unit:** 1 pair

### Indications

- Heel Spur,
- shock absorption

**Sizes:** Comes in three sizes: Small, Medium and Large



Reference	Size	Shoe size (US)	Color
GB-930070S	Small	W up to 7½ / M up to 6	Translucent/blue spot
GB-930070M	Medium	W 8 -10 / M 7-8½	Translucent/blue spot
GB-930070L	Large	W 10½-13½ - M 9-12	Translucent/blue spot


## Even Thread Rubber Wedge

Glue in place

**Packing unit:** 1 pair

**Sizes:** Comes in three sizes

Reference	Size
OD-055301	1
OD-055302	2
OD-055303	3



## Flexa Splayfoot Support

With grained sheepskin cover. Vegetable tanned. Moulded part of foam rubber. Distinct splayfoot pad. Anti-slip

**Packing unit:** 1 pair

**Sizes:** 35 to 46



**Ref.:** OD-2510XX (\*)

(\*) To order replace XX in the reference with the desired size. Example: 082510-37 for size 37

## Bio Flatfoot and splayfoot support

With grained sheepskin cover. Vegetable tanned. Super light. Flexible. Good, effective support..

**Packing unit:** 1 pair

**Sizes:** 35 to 46



**Ref.:** OD-2531XX (\*\*)

(\*\*) To order replace XX in the reference with the desired size. Example: 082531-41 for size 41

## Silicone Heel Cushion

Heel Cushion made of silicone.  
Reduces shock load onto ankle, knee and hip joints as well as onto spine.  
With raised edges.

### Indications

Shock-absorption, damping.

**Height:** 6 mm

**Packing unit:** 1 pair



Reference	Size	Shoe size (US)	Color
GB-930050S	Small	W up to 7½ / M up to 6	Translucent
GB-930050M	Medium	W 8-10 / M 7-8½	Translucent
GB-930050L	Large	W 10½-13½ - M 9-12	Translucent

## Pedisoft® Heel Cushion

### PEDISOFT® HEEL CUSHION WITH LATERAL SOFT SPOT

Viscoelastic heel insert for shock absorption and pressure relief as well as heel spur.

Packing unit: 1 pair.

### Indications

Correction, pressure protection,  
prophylaxis in case of diabetic foot.

### Special product features

The silicone material offers decisive advantages:

It is soft and smooth just like the skin and comfortable to wear.

Thanks to the memory effect it adjusts itself optimally to the contour of each part of the foot without slipping and returns to its original shape immediately again after wearing.

After use it is washable with soap and can be disinfected if required, and therefore is more hygienic than other materials.

It is long-lasting, wear-resistant and designed for multiple use  
Available in two sizes



Reference	Size	Shoe size (US)	Color
GB-107070S	Small	Women up 9, Men up to 7½	White/blue spot
GB-107070L	Large	Women Over 9, Men Over 7½	White/blue spot

## FOOT ORTHOTICS

### ORTHO KID

Purple and Black

1 PAIR PER PACK	
Reference	Size
ER-0092-25/26	25-26
ER-0092-27/28	27-28
ER-0092-29/30	29-30
ER-0092-31/32	31-32
ER-0092-33/34	33-34



### SHELL ORTHOTIC CLEAR LADIES

1 PAIR PER PACK		
Size	Reference	
	Wide	Narrow
34.5 - 35.5	AC-6403345	AC-6402345
36 - 37	AC-6403360	AC-6402360
37.5	AC-6403375	AC-6402375
38	AC-6403380	AC-6402380
39	AC-6403390	AC-6402390
39.5	AC-6403395	AC-6402395
40	AC-6403400	AC-6402400
41	AC-6403410	AC-6402410
41.5 - 42	AC-6403415	AC-6402415
42.5 - 43	AC-6403425	AC-6402425



### SHELL ORTHOTIC CLEAR MEN



1 PAIRS PER PACK	
Reference	Size
AC-6405380	38 - 39
AC-6405395	39.5 - 40
AC-6405410	41 - 41.5
AC-6405420	42 - 42.5
AC-6405430	43 - 44
AC-6405445	44.5 - 45
AC-6405455	45.5 - 46
AC-6405470	47 - 48

## FOOT ORTHOTICS

### SHELL ORTHOTIC FIRM LADIES

Color: Blue

1 PAIR PER PACK	
Reference	Size
AC-6404345	34.5 - 35.5
AC-6404360	36-37
AC-6404375	37.5
AC-6404380	38
AC-6404390	39
AC-6404395	39.5
AC-6404340	40
AC-6404341	41
AC-6404415	41.5 - 42
AC-6404425	42.5 - 43.5



### SHELL ORTHOTIC FIRM MEN

Color: Blue

1 PAIR PER PACK	
Reference	Size
AC-6406380	38 - 39
AC-6406395	39.5 - 40
AC-6406410	41 - 41.5
AC-6406420	42 - 42.5
AC-6406430	43 - 44
AC-6406445	44.5 - 45
AC-6406455	45.5 - 46
AC-6406470	47 - 48

### J-PLAS PRE-FABES



### SHELL ORTHOTIC ELEG WHITE COLOR



1 PAIR PER PACK	
Reference	Size
AC-640135	35-36
AC-640137	37-38
AC-640139	39-40
AC-640141	41-42
AC-640143	43-44
AC-640145	45-46

1 PAIR PER PACK	
Reference	Description
JM-515M07-08	Men, Size 7-8
JM-515M09-10	Men, Size 9-10
JM-515M11-12	Men, Size 11-12
JM-515W05-06	Women, Size 5-6
JM-515M07-08	Women, Size 7-8
JM-515M07-08	Women, Size 9-10